



DENSITY TESTER

DENSITY TESTER

Density Tester LDMT-A10

Labtron LDMT-A10 density tester is simple and compact. The control box adopts humanized design and the control switch adopts light-touch mode. The instrument is designed and made as per ASTM D1298 "Test Methods for Density of Petroleum and Liquid Petroleum Products (Densimeter Methods)".

Features

- Fast control rate and small over-shoot
- Adopts intelligent LCD temperature controller
- Unvarying vessel temperature
- Easy to operate and observe the sample
- Refrigeration device (Optional)

Applications

It is used to determine the density of crude petroleum and liquid petroleum.

Specifications

Model No.	LDMT-A10
Cubage of bath	Φ300*380mm
Temperature controlling range	0°C to 100°C
Temperature controlling accuracy	±0.2°C
Temperature sensor	Industrial platinum resistance, Pt100
Thermometer	Total immersion mercury thermometer
Scale division	0.2°C
Heating power	700W; 1000W
Power supply	AC220V±10%, 50Hz



Labtron Equipment Ltd.
Sentinel House, Ancells Business Park,
Harvest Crescent, Fleet, GU51 2UZ,
United Kingdom
Tel: 01252 413773
www.labtron.com • info@Labtron.com